Topmate

ELECTRIC WHEELCHAIR

EW11

<http://topmate.cc/user-manual/EW11.pdf>

<http://topmate.cc/user-manual/EW11.pdf>

<http://topmate.cc/user-manual/EW11.pdf>

<http://topmate.cc/user-manual/EW11.pdf>

CATALOGUE

1 . Product Introduction	1
2 . Assembly Instructions	2
3 . Folding Instructions	3
4 . Technical Parameter	3
5 . Joystick controller operation	4
6 . Regarding the Battery	7
7 . Pre-Driving Preparation Exercises	9
8 . Electric Wheelchair Warranty Terms	11

Welcome aboard your new power wheelchair.

This model has been designed with your practical needs in mind, equipped with modern high-tech electronics and special features for the most comfortable ride.

This manual contains important information regarding the safe use of the power wheelchair.

Please read it carefully and ensure you understand all the instructions before attempting to use the vehicle.

1.Product Introduction

Foldable armrests

Push Handles

Joystick controller

Battery

Motor

Rear wheel

Footrest

Front wheel

Scan to open YouTube and watch a tutorial on wheelchair assembly.

<https://youtu.be/BfrUekxQjvc>

2.Assembly Instructions

Note: Please care your hands safety when take folding or unfolding operation.


1.Installing the Controller Mount the controller onto the controller connector.

2.Adjust the orientation of the controller, then tighten the securing bolt.

3.Connecting the Power (The controller can be installed on the right or left side of the handlebars according to the driver's preference)

Push the wheelchair seat down to unfold it, then press the backrest down firmly to lock it in place and secure the wheelchair.

Pull the anti-tip wheel lock pin, adjust the anti-tip wheel to the correct position (inserted or extended), then release the pin, the wheels will automatically block.

Connect the power plug of the battery located under the seat.

D: Electric Mode (Powered operation)

N: Manual Mode (Freewheel operation)

Ensure that both clutch levers are in the same position before use.

3.Folding Instructions

Release the folding lock latch.

Lift the electric wheelchair seat upwards to complete the folding process.

Scan to open YouTube and watch a tutorial on wheelchair unfolding and folding.

<https://youtu.be/q05bvMtHfgI>

4.Technical Parameter

Model	EW11		
Material	Aluminium Alloy		
Unfolded Dimension	960*580*910mm/37.8*22.8*35.8in		
Folded Dimension	375*580*750mm/14.8*22.8*29.5in		
N.W.(battery not included)	18.2Kg/40.1lb	Battery Weight	2.3Kg/5.1lb
Front Wheel	7" Honeycomb tire		
Rear Wheel	11" Honeycomb tire		
Anti-tip Wheel	2" Soundless wheels		
Loading Capacity	100kgs/220.4lb		
Motor Output Power(Max)	24V 200W Dual Brushless Motors		
Battery	24V12.5Ah(Lithium)		
Range	18km/11mile		
Forward Speed	0-7km/h/0-4mph		
Charger	AC100-240V50/60HZ Output: 24V-2A		
Turning Radius	≤900mm/≤35.4in		
Brake Type	Electromagnetic Brake		
Climbing Capability	12°		

5. Joystick controller operation

The operation of the controller


Try to avoid knocking the controller, especially the joystick. Do not beat the controller and the joystick while traveling. Do not hit the controller. When operating the wheelchair, ensure that the control system is safe and reliable.

Operating environment

Do not keep the control system and its components in a humid environment for long periods of time. If food or drink is spilled on the control system, it should be cleaned immediately.

Cleaning

Use a soft cloth to gently wipe the control system and the joystick. Do not use abrasives and similar cleaners.

On/off button

The on / off button provides power to the electronic components of the control system to provide power to the wheelchair motor. In an emergency, you can use the on / off button to stop the wheelchair.

- Turn on the power: Press the power button, Indicator will display the current battery power level, and speed.
- Turn off the power: Press the power button, All indicator light will turn off.

Battery indicator

The battery indicator shows the current battery power of the wheelchair. The indicator light tells you that the wheelchair is on. At the same time, the power display also shows the status of the wheelchair.

Horn

Press the horn to sound an alert.

Acceleration button

This button can increase the speed setting gear and the corresponding speed indicator lights on.

Speed indicator (S1-S5)

This display indicates the wheelchair setting. With a total of five settings, S1 is the lowest speed, and S5 is the highest.

Deceleration button

This button reduces the speed setting gear and turns the corresponding speed indicator lights off.

Joystick

The main function of the joystick is to control the speed and direction of the wheelchair. When you release the joystick, the wheelchair automatically brakes.

Control system status indication

The battery indicator and speed indicator shows the status of the control system. See "Controller status indication description" for details.

A: Power display is stable

The current remaining battery power is displayed according to the battery voltage.

B: V1-V5 flashes in sequence

The wheelchair is charging. Unplug the charger before you can start the wheelchair

C: The control system is turned on, accompanied by sound & light alarms, and the control system is blocked from moving the wheelchair.

This indicates that there is a problem with the system. There is a problem in the electronic system of the wheelchair. Please follow the following procedures.

- Turn off;
- Make sure all connections are in place and in good contact;
- Check the battery status;
- Turn on the control system again and try to drive. If the safety circuit works again, turn it off and do not use the wheelchair again. Follow the "self-repair instructions" for operation.

Self-help guidance

If a system failure occurs, check the position of the flashing lamp on the speed display and battery indicator. Below are a series of errors which could be indicated by an abnormal sound:

Model EW11 Controller status indication

Error indication	Fault speculation	Solution
The 1st speed indicator and 2nd battery indicator flash at same time	Hall fault (left / right motors).	Check whether the motor and controller are wired correctly.
	Electromagnetic brake- cable detached (left / right motors)	
The 1st speed indicator and the 1st, 2nd battery indicator flash at same time	Motor connector unplugged (left / right motors)	Check whether the motor and controller are wired correctly.
The 3rd speed indicator and 1st battery indicator flash at same time	Joystick not in position	Turn off the power of the controller, gently toggle the joystick to make sure there is no jamming ,then turn on the controller on again.
The 3rd speed indicator and 2nd battery indicator flash at same time	The controller joystick dropped off	Contact seller for a new joystick controller.
The 4th speed indicator and 1st battery indicator flash at same time	Electric current abnormality	Contact seller for a new joystick controller.
The 4th speed indicator and 2nd battery indicator flash at same time	Driver tube breakdown	Contact seller for a new joystick controller.
The 4th speed indicator and 3rd battery indicator flash at same time	Control board short circuit	Contact seller for a new joystick controller.
The 1st speed indicator and 5th battery indicator flash at same time, and with a beep.	Upper/lower controller communication lost	Contact seller for a new joystick controller.
The 4th speed indicator and 5th battery indicator flash at same time	Electro-magnetic drive off	Contact seller for a new joystick controller.
The 5th speed indicator and 1st battery indicator flash at same time, but no' beep.	Low battery voltage	Low battery power, please charge the battery.

<p>The 5th speed indicator flash, and 1-5th battery indicator flash at same time.</p>	<p>High battery voltage</p>	<p>Battery charge is too high (and the voltage is higher than the national standard $31.5 \pm 0.5V$), need to discharge the battery.</p>
---	-----------------------------	---

6.Regarding the Battery

How to remove the battery?

Disconnect the connector plug from the battery.

With one hand, lift the battery lock mechanism upwards, and with the other hand, grab the battery handle and pull the battery outwards.

Scan to open YouTube and watch a tutorial on wheelchair battery removal.

https://youtu.be/bU21-_HdYJs

Two charging methods

First kind: Remove the battery for charging.

Two kind: Charge the controller.

Scan to open YouTube and watch a tutorial on two wheelchair charging methods.

<https://youtu.be/N2KNQao3e-Q>

Battery use

The electric wheelchair uses a 24 voltage.

It is a good idea to carry a replacement 24V/12.5Ah battery, in case you need to charge the one in use. This is a maintenance-free sealed lithium battery.

The batteries are designed to handle deep discharges. Although they are similar in performance to car batteries, they cannot be replaced. Car batteries are not safe for use in electric wheelchairs.

Battery charge and discharge

When the indicator light of the charger turns green, unplug the plug of the charger, then pull out the power plug.

In order to maximize the efficiency of your wheelchair battery, you should:

1. Fully charge your battery during initial use.
2. Start slowly when using the wheelchair in your home or on the driveway.
3. Give your battery charge 4-6 hours, and then re-operate the wheelchair.
4. After 4-5 cycles of charging and discharging, the battery effect can reach 100%, and can continue for a long time.

Important information about the battery

Frequently charged batteries provide reliable performance and long battery life. If the battery is not charged often or not fully charged, it may cause damage to the battery permanently, causing unreliable operation and reduced battery life.

If you do not often use your electric wheelchair, you must charge once a month to maintain battery life.

If you want to ship your wheelchair, you must contact the shipping company in advance to ensure that the wheelchair meets the special requirements of the transport, such as aircraft, buses and trains.

We recommend that you contact the public transport department for packaging and shipping according to special performance.

Used batteries should be carefully discarded according to your local council's regulations.

Charger

The charger is a very important part of your wheelchair and provides energy to your wheelchair quickly and safely.

Your electric wheelchair is designed to be maintenance free, however we recommend the following checks be completed on a regular basis:

Joystick/controller:

Ensure that the controller is kept away from the damp air as this can easily damage it and void the warranty.

Battery cable:

Check the battery cables to make sure they are not corroded.

Clean your wheelchair with a damp cloth:

Do not use water to clean your wheelchair.

7.Pre–Driving Preparation Exercises

Failure to follow the instructions below may result in serious injury or damage to the electric wheelchair.

Practice before use

Before venturing onto the roads, find an open area such as a park and ask a friend to accompany you until you have mastered the power chair's operation.

The power button should always be off before you get on or off the seat.

Set the speed control button appropriate for the vicinity in which you are travelling.

We recommend keeping to a slow speed until you are familiar with handling the wheelchair.

In the event of one of the following errors, immediately turn off the electric wheelchair's power supply:

- 1.Sudden movement
- 2.Operation direction is not controlled
- 3.Brake failure

Safety instructions

- Do not carry passengers.
- Do not drive across a slope.
- Do not drive under the influence of alcohol.
- Do not tailgate another vehicle.
- Do not operate hand-held communication devices such as radios or mobile phones.
- Do not exceed the safe climbing maximum angle.
- Do not use the joystick in an erratic manner when going up or downhill.
- Do not exceed the maximum user weight.

- Do not press the power button of the controller whilst moving as this could damage the electromagnetic brakes or the joystick controller.
- Do not travel on soft terrain.
- Do not attempt to climb stairs over two inches tall unless a ramp is installed.
- Do not disassemble the power wheelchair unless the electromagnetic brakes are engaged and the joystick controller is turned off.
- Do not operate the power wheelchair without a driver.
- Do not drive the power wheelchair on the main roads – the footpath must be used.
- Do not change driving direction or reverse before the power wheelchair has completely stopped.
- Always engage a slow speed when going down gradients (move the joystick to the central position)
- Always use the safety belt.
- Always keep feet on the footrest while travelling.
- Always make sure the batteries are fully charged before setting out on a journey.
- Always charge the power wheelchair in a well-ventilated area.
- Always replace a loose screw from the centre of the wheel with a wrench.

Reminder! When using an electric wheelchair, the driver must follow all the rules and regulations of the road.

Electric Wheelchair Warranty Terms

Thank you for purchasing our electric wheelchair. To protect your rights, please read the following warranty terms carefully:

Whole Vehicle Warranty Terms

1.Warranty Period:The warranty period for the entire vehicle is 2 years from the date of purchase.

2.Warranty Coverage

- All original parts of the electric wheelchair, including but not limited to the frame, wheels, braking system, handlebars, etc.
- Faults or damages not caused by human factors.

3.Exclusions from Warranty

- Damages caused by user disassembly or modifications.
- Damages resulting from accidents, collisions, misuse, or incorrect use.
- Normal wear and tear (such as wear on tires and brake pads).

4.Warranty Service Process

- To apply for warranty service, please provide proof of purchase and the product serial number.
- Please bring the electric wheelchair to a designated service center or contact our customer service center for assistance.

Battery Warranty Terms

1.Warranty Period:The warranty period for the battery is 1 year from the date of purchase.

2.Warranty Coverage

- Battery faults or performance degradation due to manufacturing defects or material issues.
- Abnormal conditions under normal usage.

3.Exclusions from Warranty

- Damages caused by improper use or failure to follow the charging instructions in the user manual.
- Damages resulting from external physical damage (such as collisions or compression).
- Battery issues that arise after the warranty period has expired.

4.Warranty Service Process

- To apply for battery warranty service, please provide proof of purchase.
- Please bring the battery to a designated service center or contact our customer service center for assistance.

Thank you for your trust and support in our products. If you have any questions or need further assistance, please feel free to contact our customer service team

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