# Mobility Scooter Manual

## ES41

<http://topmate.cc/user-manual/ES41.pdf>

<http://topmate.cc/user-manual/ES41.pdf>

<http://topmate.cc/user-manual/ES41.pdf>

<http://topmate.cc/user-manual/ES41.pdf>

# CATALOGUE

1 . Product Introduction	1
2 . Assembly tutorial	2
3 . Battery removal & charging	3
4 . Electric mode & Manual push mode	4
5 . Folding and unfolding	5
6 . Parameters	6
7 . Product standard configuration	6
8 . Inspection and maintenance before driving	7
9 . Battery instruction	7
10 . Attentions for safe driving	8
11 . Eliminate of frequent malfunction	9
12. Friendly reminder	9
13 . Electric Scooter Warranty Terms	10

Thank you for your patronage !

- Proper use of this product is the prerequisite to ensure safety, convenience and durability
- Please read this instruction carefully before using, and be sure to master the matters marked with the mark

# 1.Product Introduction

## 2. Assembly Tutorial

1

Take the scooter and the required mounting accessories out of the package

3

Hold the steering wheel angle adjuster and adjust to the desired angle

2

Unfasten the Velcro strap

4

Pull up the seat frame

5

Unfold the backrest, then place the seat cushion on the seat frame

6

Adjust the armrest mounting seat to the appropriate position, then tighten the long handle screw. (Repeat this operation on the other side)

7

Adjust the armrest to the desired height, then tighten the long handle screw. (Repeat this operation on the other side)

8

Hold the adjustment button and adjust the headrest to the desired height

Scan to open YouTube to watch assembly tutorial

<https://youtu.be/37YFFITePVE>

## Battery Disassembly:

1. Hang the basket on the basket hook.

2. Press the basket clasp downward and unfold the basket.

3. Place the bottom of the basket flat downwards.

Scan to open YouTube to watch the tutorial on installing the basket on the scooter  
<https://youtu.be/4gP7r9qS1Y4>

## 3. Battery removal & charging

Note: Please make sure the power switch turned OFF when make battery Install/Remove operation.

### How to remove the battery?

Grasp the battery handle and lift the battery upwards to remove it.

### How to charge the scooter?


On the left side of the front of the scooter, there is a charging port. Use the standard charger to charge the scooter.

Scan to open YouTube to watch the tutorial on removing the battery and charging it  
<https://youtu.be/FGMyNvBCLd0>

## 4. Electric mode & Manual push mode


The brake handle is located at the left rear of the scooter. When the brake handle is in the down position (manual/push mode), the scooter is in a brake-free state, and it cannot be powered by the battery, only pushed. You must push the brake handle up to the up position for the scooter to be in electric mode and powered by the battery.

Please note that if you push the scooter too fast while the wheels are moving freely, the scooter will slow down.

## 5. Folding and unfolding

### Unfolding Tutorial:

### Folding instructions:

Scan to open YouTube to watch the tutorial on how to unfold and fold the scooter  
<https://youtu.be/IFUDCQK3Jqw>

## 6.Parameters

Material	Magnesium alloy
Dimension	106*45.5*90cm/41.7*17.9*35.4in
Folded Dimension	106*45.5*43.5cm/41.7*17.9*17.1in
N.W.(w/o battery)	17.7Kg/39Lb
Battery Weight	1.7Kg/3.7Lb
Loading Capacity	125Kg/275.6Lb
Motor Rated Power	24V180W
Motor Rated Speed	4000±100 rpm/min
Motor	Brushless Motor
Maximum Torque	≥25N.m
Battery	24V 10Ah 240Wh
Max Speed	0-7KM/h (0-4MPH)
Range	12-15KM/7-9MPH
Reverse Speed	3KM/h(2MPH)
Drive system	Rear-wheel drive
Charger	29.4V 2A
Charging time	5-6H
Brake Type	Electromagnetic brake
Climbing Capability	0-12°

The above data are for reference only. The technical data of this product series are subject to change without prior notice. In actual use, due to different road conditions, climate and driving speed, data such as charging time and driving range are different.

## 7.Product standard configuration

Portable electric folding scooter	1 unit	Charger	1 unit
Installation tool	1 unit	Instructions	1 unit
Car Basket	1 unit		

## 8. Inspection and maintenance before driving

In order to ensure that this portable electric folding car can bring convenience to your life, please check and maintain the product according to the following requirements before using.

1. Check whether the battery is full of electricity to ensure that all kinds of cables are properly connected.
2. Check whether the brake system is sensitive and effective.
3. Check whether the tire pressure is normal
4. Before driving, please correct the front wheel and tighten the handle fastening handle.
5. The handlebar and the moving clearance of each part, should be inspect fastening and lubricate irregularly

## 9. Battery Instruction

★ Improper use will shorten the service life of the battery, and may lead to the fever of the battery, Fire, or even explosions. Please read the following notes carefully:

★ Do not separate and decompose the battery. There are protective circuit function parts inside the battery pack. Disassembly will damage these parts;

★ It is strictly prohibited to short circuit the positive and negative terminals of the battery, Once a short circuit is found, it will cause the battery heat, fire and other dangers;

★ Do not burn or heat the battery. Heating or burning the battery will melt the insulation, thus causing the battery fire, explosion and other dangers;

★ Do not immerse the battery in water or in a humidity environment, The water in the battery pack is easy to cause the short circuit of the protective circuit board in the battery pack, Thus causing short circuit heating and other dangers;

★ You must charge the battery with a charger that matches the battery, Otherwise, abnormal chemical reactions will occur in the battery, causing the battery to heat up, burn and other dangers;

★ It is strictly prohibited to use a hammer or other hard objects to hit the battery, cutting the battery with sharp tools, When the battery is damaged, it may cause short circuit, thus causing fire and other dangers;

★ If the battery is damaged or leaked, and the leakage material enters the eyes carelessly, please don't rub it with your hands, wash it with clean water immediately, and go to the doctor as soon as possible;

★ Do not mix lithium ion batteries made by different manufacturers or with different capacity. If so, it may cause heat, fire and other dangers in the process of charging and discharging;

★ The battery needs to be fully charged in time (within 3 days) after use, If the battery pack is not used for a long time, it must be fully charged every two months, otherwise the battery will be damaged easily;

- ★Do not exceed the battery charging time (generally within 10 hours), for fear that the danger of overheating. If the battery charging exceeds the specified time and is not fully charged, please stop charging;
- ★Don't put the battery where a child can touch so as to avoid danger;
- ★Do not store the battery in microwave or high pressure environment so as to avoid danger.

## 10. Attention for safe driving

In order to ensure that you can use this portable electric folding car safely, please be sure to follow the following requirements for proper operation:

- ★Before driving for the first time, be sure to do driving exercises in the open field. Only drive on the road after mastering acceleration, deceleration, braking, honking and other performances perfectly
- ★Follow traffic rules and maintain moderate speed; Driving at top speed will not only dramatically increase power consumption and reduce range, but also reduce safety
- ★When driving in a crowded indoor environment, such as a supermarket, you must keep a low speed to control the braking distance. When turning, the speed must not exceed 2km/h
- ★It should not be used in rainy days, and wading should be avoided as far as possible. When wading through, the wading depth should not exceed 20mm and the speed should not exceed 1km/h
- ★With the power switch key, in the non-driving state, the power switch must be closed and the key pulled out, to prevent children and others from wrong operation

## 11. Eliminate of frequent malfunction

<b>1 Beep</b>	The Battery needs to be charged or there is a bad connection. Check the battery connections, if the battery connections are good, charge the battery.
<b>2 Beep</b>	There is a bad connection to the motor. Check all connection between the motor and the controller
<b>3 Beep</b>	The motor has a short circuit of a battery connection. Please contact your service agent.
<b>4 Beep</b>	Not in use.
<b>5 Beep</b>	Not in use.
<b>6 Beep</b>	The charging corf of the charging port is not pulled out.
<b>7 Beep</b>	Controllever is out of order. Please check whether the control lever is in the normal position; Make sure not to touch the controllever before turning on the power.
<b>8 Beep</b>	A controller fault is indicated. Make sure al connections are secure.
<b>9 Beep</b>	The position of the manual brake handle is incorrect. Please put the brake handle upward.
<b>10 Beep</b>	Battery voltage is too high and need to call the service agent.

Please note that all sounds are based on the scooter's buzzer, if the above instructions do not work, please contact your agent for judgment and repair to avoid secondary damage to the scooter!

## 12. Friendly reminder

★When the scooter starts, please slow down the throttle to prevent the front wheel from skidding

★When the scooter climbing the slope, please lean forward to Keep the center of gravity steady and prevent scooter from turning over after climbing the slope

★When the scooter is turning, please slow down and avoid turning over

★in the normal using of the scooter, please check whether the screw loose or off periodically

★Please drive as little as possible when the road is slippery in rainy days

# Electric Scooter Warranty Terms

Thank you for purchasing our electric scooter. To protect your rights, please read the following warranty terms carefully:

## Whole Vehicle Warranty Terms

1. Warranty Period: The warranty period for the entire vehicle is 2 years from the date of purchase.

### 2. Warranty Coverage

- All original parts of the electric scooter, including but not limited to the frame, wheels, braking system, handlebars, etc.
- Faults or damages not caused by human factors.

### 3. Exclusions from Warranty

- Damages caused by user disassembly or modifications.
- Damages resulting from accidents, collisions, misuse, or incorrect use.
- Normal wear and tear (such as wear on tires and brake pads).

### 4. Warranty Service Process

- To apply for warranty service, please provide proof of purchase and the product serial number.
- Please bring the electric scooter to a designated service center or contact our customer service center for assistance.

## Battery Warranty Terms

1. Warranty Period: The warranty period for the battery is 1 year from the date of purchase.

### 2. Warranty Coverage

- Battery faults or performance degradation due to manufacturing defects or material issues.
- Abnormal conditions under normal usage.

### 3. Exclusions from Warranty

- Damages caused by improper use or failure to follow the charging instructions in the user manual.
- Damages resulting from external physical damage (such as collisions or compression).
- Battery issues that arise after the warranty period has expired.

### 4. Warranty Service Process

- To apply for battery warranty service, please provide proof of purchase.
- Please bring the battery to a designated service center or contact our customer service center for assistance.

Thank you for your trust and support in our products. If you have any questions or need further assistance, please feel free to contact our customer service team

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